

# Setting The New Standard



## The MK1 Manual Keratometer

# The Reichert MK1 Keratometer... an industry standard

At Reichert, we continuously strive to innovate and improve all our products. So how do you improve an instrument that has been the standard for more than half a century? Take away the need to replace bulbs by using light emitting diodes (LEDs).

The Reichert MK1 Keratometer now utilizes LED technology to provide clear, consistent lighting, removing the need to ever change bulbs. What we didn't change are all the great features that have made our Keratometer the worldwide instrument of choice by discerning practitioners for over fifty years:

- Both major and minor meridians can be measured once major axis is located
- Accurate measurement of corneal curvatures
- Allows subjective review of corneal surface
- Measures concave surfaces for contact lens fitting
- Requires no patient response or judgement
- Easy elevation and lowering allows quick patient positioning
- Coincidence focusing system
- Engraved focusing scale on eyepiece
- Dual eye level sighting for vertical measurement
- Single position measurement in both meridians
- Dial readings are shown in both diopters and millimeters for maximum convenience, covering both corneal and contact lens measurement
- Clear, uniform illuminated mires
- Accurate to .012D from 36D to 52D
- Engraved focusing scale on eyepiece

**Call your authorized Reichert Ophthalmic Instruments Distributor today.**

Distributor list available at [www.reichertoi.com](http://www.reichertoi.com)

## Reichert Quality

The Reichert name ensures a high quality, reliable instrument that you can count on every time.



## Catalog Numbers & Specifications

<b>15040</b>	MK1 Manual Keratometer
<b>Power Range</b>	36D to 52D
<b>Radius Range</b>	6.50mm to 9.40mm
<b>Accuracy</b>	0.12D
<b>Height</b>	15.5 in (39.4 cm)
<b>Length</b>	18.5 in (39.4 cm)
<b>Width</b>	6.75 in (17.1 cm)

*Specifications subject to change without notice*

## Important Notice

For your own protection, only purchase your Reichert Ophthalmic Instruments from, and have them serviced by, an Authorized Reichert Distributor.



**Reichert**  
Ophthalmic Instruments

Reichert, Inc.

3374 Walden Avenue, Depew, NY 14043  
Toll Free: 888-849-8955 | Tel: (716) 686-4500  
Fax: (716) 686-4545 | [www.reichert.com](http://www.reichert.com)